

Product Datasheet

Product code: 3.DX.0

Product type: Biology Tweezers

Made in: Switzerland

Description: Type 3 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The high Alloy Austenitic Stainless Steel material provides improved performance in severe environments.

Tip description: straight, ultra fine, pointed, superior finish

Coating/finishing: polished

Body material: Anti-Magnetic Anti-Acid High-alloy Stainless Steel (DX)


Picture:




Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	12.0 mm – 0.472 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	0.08 mm – 0.003 inch
Tip thickness	0.04 mm – 0.002 inch