

# Product Datasheet

**Product code:** **71MZ.SA.0**

**Product type:** Ceramic Replaceable Tip Tweezers

**Made in:** Switzerland

**Description:** Type 71MZ Ceramic replaceable tip tweezers are suitable to applications which require holding a small component firmly. The Ceramic tips provide the tweezers with outstanding resistance to high temperatures and aggressive environments combined with exceptional hardness.

**Tip description:** straight, fine, beveled, superior finish

**Tip material:** Zirconia Toughened Alumina (MZ)

**Coating/finishing:** Sandblasted


**Body material:** Anti-Magnetic Anti-Acid Stainless Steel (SA)

**Picture:**



**Tweezer dimensions:**

OAL	135.0 mm – 5.315 inch
Width	9.0 mm – 0.354 inch
Height	12.0 mm – 0.472 inch



**Tip dimensions:**

Tip width	0.70 mm – 0.028 inch
Tip thickness	0.60 mm – 0.024 inch