

Product Datasheet

Product code: **MPTSS3**

Product type: Stainless steel spatulas


Made in: Pakistan

Description: MPTSS3 Stainless Steel spatula is suitable to collect material. With its double tip rounded body, the MPTSS3 is an auxiliary tool for laboratory use. The spoon tips consist in a short, flat, squared tip and a long, flat, rounded tip. Stainless Steel probes are rust-free and equipped with knurled handles which prevent any possible slips.

Tip description: short and flat with squared tip, long flat with rounded tip

Body material: Stainless Steel (S)

Picture:



Product dimensions:

OAL	165.0 mm – 6.496 inch
Width	9.5 mm – 0.374 inch
Height	6.0 mm – 0.236 inch

Tip (1) details:



Tip (1) dimensions:

Tip width	8 mm – 0.315 inch
Tip thickness	0.5 mm – 0.020 inch

Tip (2) details:



Tip (2) dimensions:

Tip width	8.7 mm – 0.343 inch
Tip thickness	0.5 mm – 0.020 inch