

## Product Datasheet

**Product code:** **2WL.SA.1**

**Product type:** Wafer Tweezers

**Made in:** Switzerland

**Description:** Type 2WL wafer tweezers are suitable to delicate applications which require holding a 2-inch flat component firmly. The line serration provides the handles with extra grip. The bent paddle helps to hold the component still during positioning.

**Tip description:** straight, wafer

**Coating/finishing:** Sandblasted

**Body material:** Anti-Magnetic Anti-Acid Stainless Steel (SA)


**Picture:**




**Tweezer dimensions:**

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	12.0 mm – 0.472 inch

**Tip details:**



Dimensions inch



Dimensions mm

**Tip dimensions:**

Tip width	6.20 mm – 0.244 inch
B	9.00 mm – 0.354 inch
C	2.40 mm – 0.094 inch
D	3.30 mm – 0.130 inch
E	2.30 mm – 0.091 inch