

Product Datasheet

Product code: **5.SA.GP.0**

Product type: Gold plated tweezers

Made in: Switzerland

Description: Type 5 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The coating of 2-microns thickness, made of pure 24-carat gold plating, provides excellent scratch resistance and makes them particularly suitable for handling nickel grids in immuno-gold work. The tweezers exhibit remarkable durability under high-volume autoclaving conditions, and the gold plating enhances their resistance to chemical corrosion, oxidation, and electrical conductivity.

Tip description: straight, very fine, sharp, superior finish

Coating/finishing: Gold plated (GP)

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)


Picture:




Tweezer dimensions:

OAL	110.0 mm – 4.331 inch
Width	10.0 mm – 0.394 inch
Height	12.0 mm – 0.472 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	0.10 mm – 0.004 inch
Tip thickness	0.06 mm – 0.002 inch