

Product Datasheet

Product code: **85C.SA.1**

Product type: Wafer Tweezers

Made in: Switzerland

Description: Type 85C wafer tweezers are suitable to delicate applications which require holding a 2-inch flat component firmly. The line serration provides the handles with extra grip. The bent paddle helps to hold the component still during positioning.

Tip description: straight, wafer

Coating/finishing: Sandblasted

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)


Picture:




Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	12.0 mm – 0.472 inch

Tip details:



Dimensions inch



Dimensions mm

Tip dimensions:

Tip width	6.50 mm – 0.256 inch
B	4.50 mm – 0.177 inch
C	2.90 mm – 0.114 inch
D	3.90 mm – 0.154 inch
E	2.00 mm – 0.079 inch