

Product Datasheet

Product code: **MPTSS1**

Product type: Stainless steel spatulas


Made in: Pakistan

Description: MPTSS1 Stainless Steel spatula is suitable to collect material. With its double tip rounded body, the MPTSS1 is an auxiliary tool for laboratory use. The spoon tips consist in a short, pointed tip and a long, curved, rounded tip. Stainless Steel probes are rust-free and equipped with knurled handles which prevent any possible slips.

Tip description: pointed tip arrow-shaped, flat long curved with a rounded tip

Body material: Stainless Steel (S)

Picture:



Product dimensions:

OAL	180.0 mm – 7.087 inch
Width	7.8 mm – 0.307 inch
Height	6.0 mm – 0.236 inch

Tip (1) details:



Tip (1) dimensions:

Tip width	0.6 mm – 0.024 inch
Tip thickness	0.5 mm – 0.020 inch

Tip (2) details:



Tip (2) dimensions:

Tip width	5.3 mm – 0.209 inch
Tip thickness	0.8 mm – 0.031 inch